

SAFETY DOUBLE GRAB BAR - 929 mm



| | |
|---------------|----------------|
| Product Code | 28974 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | No |
| CH-REP | Yes |

| | |
|------------|----------------|
| RDM (NSIS) | 2535067 |
| CND | Y181899 |
| EAN/UPC | 8023279289749 |
| PARAF | 925875344 |
| GTIN | 08023279289749 |
| GMDN | 35584 |

Description

SAFETY DOUBLE GRAB BAR - 929 mm

Innovative line of GRAB BAR, can be quickly installed everywhere without need of drilling and screwing.
Patented suction cups (Ø 93 mm).
Holds up to 85 kg vertical and 60 kg lateral weight (tested by SGS).

Ideal solution to support elder people during daily bathroom use.
Top quality product made of heavy duty ABS, with patented "green/red" suction indicator for better safety.

Multilingual box and manual: GB, FR, IT, ES, PT, DE, GR, Arabic.