

# FIXED PLASTIC GRAB BAR - 60 cm



Product Code	28977
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required

EAN/UPC	8034055100646
PARAF	977812864

## Description

### PLASTIC RIBBED GRAB BAR - 60 cm

Ribbed plastic grab bar in Medical PVC, complete with ABS mounting plates and screw covers.

Wall screws not included.

Electrically isolated.

The ribbed bars offer excellent anti-slip grid and are warm to the touch.

Can be mounted either horizontally or vertically.

Tube Ø 3.5 cm.

Max user weight 120 kg (265 lb).

Manual: GB, IT.

Made in Italy.