


U-SHAPE POWDER COATED STEEL FOLDABLE BAR - 70 cm



Product Code	28979
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2563436
CND	V0880
EAN/UPC	8034055101735
PARAF	977812876

Description

U-SHAPE POWDER COATED STEEL FOLDING BAR - 70 cm

White powder coated wall mounted lift up toilet support.
Designed to help getting up or transferring to and from the toilet.
It can be folded upwards for easy accessibility.

Max user weight 125 kg (275 lb).

Wall screws not included.

Manual: GB, FR, IT, ES, DE, GR.

Made in Italy.