



Product Code	49021
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

## Description

### LUXIFLEX II PLUS on trolley - LED

Powerful LED light source that offers exceptional light quality. Flexible, height-adjustable and highly precise arm allows adjustment to any position. Compact head design and switch located at the base of the lamp. Class II.

Led Plus II with CRI 97%

Plus II has the high performance of the color Ra (97%) that depicts the exceptional quality of white light for an excellent distinction of fabrics and reproduction of all visible colours, naturally and precisely.

Manual in: GB, FR, IT, DE, ES, PT

Made in Spain.

## Technical Specifications

Consumption (Watt) 12  
 Life of the bulbs (hours) 50,000  
 Light intensity at 50 cm (lux) 50,000  
 Light intensity at 1 m (lux) 17,000  
 Colour temperature  $\pm 5\%$  at 1 m (K) 4,500°  
 Colour rendering index (Ra)% 97  
 Light field size at 50 cm ( $\varnothing$  mm) 150  
 External  $\varnothing$  of the reflector (mm) 65  
 Colour white  
 Norms EN 60601-1-2  
 Marks EU UE 2017/745  
 Weight (kg): -  
 Operating voltage V 100/240  
 Operating Hz 50/60