

LED SPOT-ADJUSTED LIGHT - trolley



Product Code	49047
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
EAN/UPC	8023279490473

Description

LED SPOT-ADJUSTED LIGHT - trolley

LED Spot examination lamp with multi-optical lens and light spot adjustable by inner mini light grating.

- Additional features:
- uniform and pure light spot
 - front optical lens
 - inner light grating
 - rear dual optical lens
 - 100 cm super long 360° free position and PVC coated flexible metal arm

Manual: GB, FR, IT, ES, PT, DE, PL, GR.

Stainless steel handle.
Super long flexible soft arm, length 100 cm.
Light spot adjustable by circle on lamp head.
Digital brightness control.
ABS module base.

Technical Specifications

Bulbs/power: 1 LED / 5 W
Ø of lamp head: 25 mm
Illumination:
25,000 lux at 30 cm
10,000 lux at 50 cm
Light field Ø: max 70 mm, min 15 mm at 30 cm
Colour temperature: 4,800°-6,000° K
Lamphead: Ø 25 mm, lenght 125 mm
Power supply: 100-260 V, 50/60 Hz
Power cord: lenght 272 cm