


PREGNANCY TEST - strip 4 mm - professional



Product Code	29091
Unit of sale	box of 50 pcs
Minimum order	1
Type	IVD - Professional Use
Class	

RDM (NSIS)	2084473
CND	W0102160302
EAN/UPC	8023279290912

Description

PREGNANCY TEST 20 mIU/ml

- PROFESSIONAL TEST
- Type: strip 4 mm

GIMA offers a complete line of pregnancy tests, both for otc market (self test) and for professional use.

All products are packed in cardboard boxes with multilingual user manual (GB, IT, ES, FR, GR, PT).

Main features:

- each test is individually packed in sealed aluminium pouch
- very quick, in 3 minutes
- early detection: it can be used since the first day you miss your period
- reliability: 99.9%
- sensitivity: 20 mIU/ml
- clinically tested
- shelf life 36 months