

PH VAGINAL TEST - strip



Product Code	29110
Unit of sale	box of 10 pcs
Minimum order	1
Type	IVD - Professional Use
Class	

RDM (NSIS)	1691103
CND	W0101060299
EAN/UPC	8023279291100
PARAF	984599504

Description

Quick test to determine vaginal pH. An easy routine test to check vaginal pH, results in just 1 minute, to make sure It is within the physiological limits (< 4.5).
An increase in the pH level may indicate a change in the vaginal environment. This could mean there is a possible vaginitis. The test is ideal for pregnant women and for women with particular diseases such as diabetes.
For professional use.

Standard accessories

The kit contains: 10 cotton swabs + 10 reagent strips in an individual box with user manual (English + Italian)