


VAGINITIS TEST KIT - professional - card



Product Code	29111
Unit of sale	1 pc
Minimum order	1
Type	IVD - Professional Use
Class	

RDM (NSIS)	1691108
CND	W0101060299
EAN/UPC	8023279291117

Description

A unique, easy-to-use and revolutionary test for gynecologists and operators. Finds the causes of the inflammation/vaginal infection in just 10 minutes. It replaces traditional bacterial tests that take days and costs only a tenth. When necessary dilute vaginal fluid specimen in a buffer solution and then transfer to reagent pad. The instruction manual contains a result guide and can be used to recognize: Trichomoniasis, negative gram bacteria, Gardnerella, Chlamydia, serious bacterial vaginitis. For professional use.

WARNING: SHORT EXPIRY DATE FOR TESTS. Tests have 24 months expiry at production time. Our stock may vary from 10 to 22 months.

Standard accessories

- 1 test
- Vaginal swab
- Fluid specimen (buffer solution) in single box
- Manual (english, italian)