

OVULATION TEST - self test - strip 4mm



Product Code	29125
Unit of sale	box of 25 pcs
Minimum order	1
Type	IVD - Self Test
Class	

RDM (NSIS)	2149093
CND	W0102160304
EAN/UPC	8023279291254
PARAF	983277916

Description

OVULATION SELF TEST

- strip 4 mm
- box of 25 pcs

Easy to use self test for determination of the highest peak of fertility in woman.
Early detection test, 99,9% accurate and clinically tested.
Sensitivity: 30 mIU/ml.

All products are packed in single aluminium pouches in a cardboard boxes with multilingual user manual (GB, IT, ES, FR, GR, PT, DE, PL).

Tests have 24 months expiry at production time.
Our stock may vary from 10 to 22 months.

Results in 5 minutes.