

TRUSYNTH FAST ABSORB. SUTURE gauge 3/0 circle
3/8 needle 24mm - 70cm - undyed



Product Code	22916
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	1781318
CND	H010101010202
EAN/UPC	8023279229165
GMDN	17471

Description

TRUSYNTH FAST - braided & coated polyglactin 910 fast
Equivalent Ethicon brand: VICRYL RAPID

- Gauge: 3/0
- Metric: 2.0
- Colour: Undyed
- Circle: 3/8
- Needle profile: reverse cutting
- Needle length: 24 mm
- Suture length: 70 cm

Technical details:

- synthetic absorbable surgical suture
- braided and coated polyglactin 910
- sterilized by ethylene oxide
- shelf life of 5 years
- absorption period: 42 days

Applications:

- general soft tissue approximation

Multilanguage manual and box: GB, FR, IT, ES, PT, DE, GR.