

TRUSYNTH FAST ABSORB. SUTURE gauge 4/0 circle  
3/8 needle 19mm - 70cm - undyed



|               |                |
|---------------|----------------|
| Product Code  | 22918          |
| Unit of sale  | box of 12 pcs  |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | III            |

|        |    |
|--------|----|
| UK-REP | No |
| CH-REP | No |

|            |               |
|------------|---------------|
| RDM (NSIS) | 1781318       |
| CND        | H010101010202 |
| EAN/UPC    | 8023279229189 |
| GMDN       | 17471         |

Description

TRUSYNTH FAST - braided & coated polyglactin 910 fast  
Equivalent Ethicon brand: VICRYL RAPID

- Gauge: 4/0
- Metric: 1.5
- Colour: Undyed
- Circle: 3/8
- Needle profile: reverse cutting
- Needle length: 19 mm
- Suture length: 70 cm

Technical details:

- synthetic absorbable surgical suture
- braided and coated polyglactin 910
- sterilized by ethylene oxide
- shelf life of 5 years
- absorption period: 42 days

Applications:

- general soft tissue approximation

Multilanguage manual and box: GB, FR, IT, ES, PT, DE, GR.