

ELASTOMULL HAFT COHESIVE BANDAGE 20 m x 8 cm



Product Code	59218
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2294828
CND	M0303010201
EAN/UPC	4042809022032

Description

ELASTOMULL® COHESIVE BANDAGES
Elastic retention bandage. It provides wound care splint retention and light support.
Comfortable non woven fabric (cotton, viscose and polyamide) coated with syntethic latex.
Stretched length 20 m.

- BSN code: 45477-00
- Elastomull® Haft 20 m x 8 cm
- Box of 1 pc.
- Multilanguage box: GB, FR, IT, ES, DE, PT, NL, SE.