

ELASTOMULL HAFT COHESIVE BANDAGE 20 m x 10 cm



|               |                |
|---------------|----------------|
| Product Code  | 59219          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|            |               |
|------------|---------------|
| UK-REP     | No            |
| CH-REP     | No            |
| RDM (NSIS) | 2294821       |
| CND        | M0303010201   |
| EAN/UPC    | 4042809022063 |

Description

ELASTOMULL® COHESIVE BANDAGES  
Elastic retention bandage. It provides wound care splint retention and light support.  
Confortable non woven fabric (cotton, viscose and polyamide) coated with syntethic latex.  
Stretched length 20 m.

- BSN code: 45478-00
- Elastomull® Haft 20 m x 10 cm
- Box of 1 pc.
- Multilanguage box: GB, FR, IT, ES, DE, PT, NL, SE.