

# GIMASONIC POCKET DOPPLER WITH COLOUR DISPLAY without probe



|               |                |
|---------------|----------------|
| Product Code  | 29480          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | II A           |

|            |               |
|------------|---------------|
| RDM (NSIS) | 2681945       |
| CND        | Z12080103     |
| EAN/UPC    | 8023279294804 |
| UMDNS      | 11696         |

## Description

### COLOUR SCREEN POCKET DOPPLER

GIMASONIC Colour screen doppler - without probe

- built-in speaker
  - TFT-LCD, real-time display of FHR
  - three working modes: real-time FHR, average FHR and manual calculation
  - display of FHR, bar graph and heartbeat waveform, red alarm for abnormal FHR
  - low power consumption, continuously work for more than 8 hours with two "AA" batteries
  - indication for probe-off and low battery
  - range of probes for foetal or vascular use: 2, 3, 8 MHz probe
- Suitable for use after the 12th week of pregnancy  
Multilingual box and manual: GB, FR, IT, ES, PT, DE.

## Technical Specifications

- Ultrasound working frequency: 2 MHz  $\pm$  10%
- Fetal heart rate:
  - range: 50-240 bpm
  - resolution: 1 bpm
  - accuracy:  $\pm$  2 bpm
- Alarm:
  - high limit: 160 bpm
  - low limit: 120 bpm
- Ultrasound intensity: < 5 mW/cm<sup>2</sup>
- Auto power off after 1 minute no signal
- Dimension: 135x92x29 mm
- Weight: 245 g
- Power supply: 2 AA Batteries LR6 - 1.5 V DC