

SUTURE KIT 2 - sterile



Product Code	26949
Unit of sale	1 kit
Minimum order	1
Type	Medical device
Class	I S

RDM (NSIS)	2394559
CND	V0599
EAN/UPC	3800232424646
PARAF	963137120

Description

- Kit Content:
- Kocher forceps
 - Plastic forceps
 - Haemostatic clamp
 - Metal forceps with plastic ends
 - Littauer scissors
 - Gauze swabs 10x10 cm 8 ply
 - 4 or 5 Gauze plugs
 - WNW (TNT)+PE drape 50x60 cm
 - WNW (TNT)+PE drape 50x60 cm with hole
 - 2 adhesive WNW (Woven Non Woven) strips
 - PE bag minigrip