

SONOTRAX LITE POCKET DOPPLER without probe



|               |                |
|---------------|----------------|
| Product Code  | 29490          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | II A           |

|            |               |
|------------|---------------|
| CH-REP     | No            |
| RDM (NSIS) | 2958300       |
| CND        | Z12080103     |
| EAN/UPC    | 6944413800120 |
| UMDNS      | 11696         |
| GMDN       | 34040         |

Description

SONOTRAX LITE DOPPLER - without probe - with bag

Technical Specifications

Size: 32x85x138 mm  
Weight: 290 g (including battery)  
Heart rate range: 50-210 bpm  
Resolution: 1 bpm  
Accuracy: ± 3 bpm  
LCD display: 45x25 mm (only code 29491-93)  
Auto average and manual mode  
Recording length: 480 s / 240 s  
Auto power-off to save power  
Manual: GB, IT - on request: DE, ES, FR, PT.