

SONOTRAX PRO POCKET DOPPLER WITH DISPLAY

without probe



Product Code	29493
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

CH-REP	No
RDM (NSIS)	2640520
CND	Z12059005
EAN/UPC	6944413800175
UMDNS	11696
GMDN	34040

Description

SONOTRAX PRO II DOPPLER with backlighted display built-in audio recorder and rechargeable batteries with battery charger - without probe - with carrying bag

Technical Specifications

Size: 32x85x138 mm
Weight: 290 g (including battery)
Heart rate range: 50-210 bpm
Resolution: 1 bpm
Accuracy: ± 3 bpm
LCD display: 45x25 mm (only code 29491-93)
Auto average and manual mode
Recording length: 480 s / 240 s
Auto power-off to save power
Manual: GB, IT - on request: DE, ES, FR, PT.