

FOETAL MONITOR PROBE 1 MHz for 29506 Gima Baby Sound



Product Code	29508
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	1732835
CND	Z12080180
EAN/UPC	6941900605084

Description

1 MHz FOETAL MONITOR PROBE

The Fetal Doppler (code 29506) can be used as a portable Fetal Monitor if it works with 1MHz Fetal Monitor probe. It is suitable for long-term monitoring and will be more accurate than with the Doppler probe.