

TRUSYNTH ABSORB. SUTURE gauge 4/0 circle 1/2
needle 19mm - 75cm - violet



Product Code	22955
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III
UK-REP	No
CH-REP	No
RDM (NSIS)	1781315
CND	H010101010202
EAN/UPC	8023279229554
GMDN	17471

Description

TRUSYNTH - braided & coated polyglactin 910
Equivalent Ethicon brand: VICRYL

- Gauge: 4/0
- Metric: 1.5
- Colour: Violet
- Circle: 1/2
- Needle profile: taper cutting
- Needle length: 19 mm
- Suture length: 75 cm

- Technical details:
- braided and coated polyglactin 910
 - coating of calcium stearate and polyglactin 370
 - sterilized by ethylene oxide
 - shelf life of 5 years
 - wound support: 28-30 days
 - mass absorption: 56-70 days

- Applications:
- ophthalmic procedures
 - mucosal closures
 - skin closure
 - dental surgery
 - orthopedic surgery
 - general soft tissue approximation

Multilanguage manual and box: GB, FR, IT, ES, PT, DE, GR.