

EVOLVE DUAL HEAD STETHO - shiny head - glittered Y ceil blue



Product Code	49563
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2513128
CND	C9005
EAN/UPC	8023279495638
UMDNS	13755
PARAF	987471923
GTIN	08023279495638
GMDN	13755

Description

GIMA EVOLVE DUAL STETHOSCOPES - length 75 cm

- SHINY HEAD SERIES - glitter tubing
- Gima Evolve - ceil blue
- Binaural: mirrored
- Ring: transparent

High acoustics and sensitivity stethoscope line with 4 different heads and attractive selection of coloured tubing. Chestpiece combination with floating diaphragms which allow to hear low and high frequency sounds applying light or firm pressure on chestpiece.

Stainless steel chestpiece can catch lower frequency sound sensitively. Its inner side is polished making it smooth and reducing sound's attenuation during sound detection.

Attractive selection of fashion coloured and trendy tubes.

New electrical plating technology for longer life and corrosion resistance.

Features

- dual head chestpiece (Ø 48 mm and 36 mm) with anti-chill rings
- aluminium binaural with spring inserted into the tube (stainless steel for Rose Gold version)
- ergonomic and soft silicone eartips for best comfort
- PVC tube Ø 10 mm, length 60 cm
- total length 75 cm

Supplied with name tag and spare kit in transparent case including 2 pairs of soft eartips and 2 floating diaphragms.

Multilingual box and manual: GB, FR, IT, ES, PT, DE, NL, PL, CZ, BG, RO, HR, GR, SE, Arabic.