



Product Code	49567
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2513128
CND	C9005
EAN/UPC	8023279495676
UMDNS	13755
PARAF	987472026
GTIN	08023279495676
GMDN	13755

## Description

GIMA EVOLVE DUAL STETHOSCOPES - length 75 cm

- ALL BLACK HEAD SERIES
- Gima Evolve - black
- Binaural: black
- Ring: black

High acoustics and sensitivity stethoscope line with 4 different heads and attractive selection of coloured tubing. Chestpiece combination with floating diaphragms which allow to hear low and high frequency sounds applying light or firm pressure on chestpiece.

Stainless steel chestpiece can catch lower frequency sound sensitively. Its inner side is polished making it smooth and reducing sound's attenuation during sound detection.

Attractive selection of fashion coloured and trendy tubes.

New electrical plating technology for longer life and corrosion resistance.

### Features

- dual head chestpiece (Ø 48 mm and 36 mm) with anti-chill rings
- aluminium binaural with spring inserted into the tube (stainless steel for Rose Gold version)
- ergonomic and soft silicone eartips for best comfort
- PVC tube Ø 10 mm, length 60 cm
- total length 75 cm

Supplied with name tag and spare kit in transparent case including 2 pairs of soft eartips and 2 floating diaphragms.

Multilingual box and manual: GB, FR, IT, ES, PT, DE, NL, PL, CZ, BG, RO, HR, GR, SE, Arabic.