



Product Code	59583
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	II B

RDM (NSIS)	2477446
CND	M040299
EAN/UPC	4054596953215
PARAF	983199148

Description

KERRAMAX CARE™ SUPER-ABSORBENT DRESSING

3M™ code: PRD500-240
 Kerramax™ Care 20x22 cm
 Minimum order: Box of 10

The unique wicking layer distributes fluid both horizontally and vertically throughout the dressing, utilizing the full absorption capacity. Exu-Safe™ Technology, built into the dressing core, is an advanced method of absorbing high level of exudate from wounds, including bacteria and MMPs, creating an optimal environment for wound healing. It protects delicate wound tissue and surrounding, helping to reduce the risk of maceration, even under compression. Its heat-sealed border prevents exudate leakage.

It can be left in place for up to 7 days.
 Absorption capacity: 0.95 g/cm².
 Indicated for: pressure injuries venous leg ulcers and diabetic foot ulcers.

Latex free.
 Ethylene oxide sterilization.