

KERRAMAX CARE 3M 20x30 cm



Product Code	59584
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	II B

RDM (NSIS)	2477448
CND	M040299
EAN/UPC	4054596953307
PARAF	988754990

Description

KERRAMAX CARE™ SUPER-ABSORBENT DRESSING

3M™ code: PRD500-380-B10
Kerramax™ Care 20x30 cm
Minimum order: Box of 10

The unique wicking layer distributes fluid both horizontally and vertically throughout the dressing, utilizing the full absorption capacity. Exu-Safe™ Technology, built into the dressing core, is an advanced method of absorbing high level of exudate from wounds, including bacteria and MMPs, creating an optimal environment for wound healing. It protects delicate wound tissue and surrounding, helping to reduce the risk of maceration, even under compression. Its heat-sealed border prevents exudate leakage.

It can be left in place for up to 7 days.
Absorption capacity: 0.95 g/cm2.
Indicated for: pressure injuries venus leg ulcers and diabetic foot ulcers.

Latex free.
Ethylene oxide sterilization.