

# TEGADERM 3M SILICONE FOAM BORDER DRESSING 15x17 cm



Product Code	59598
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	II B

RDM (NSIS)	2659415
CND	M04040602
GTIN	04054596845299

## Description

### 3M™ TEGADERM™ SILICONE FOAM

3M™ code: 90647

Tegaderm™ silicone foam border dressing 15x17 cm - sacrum

3M™ Tegaderm™ Silicone Foam Border Dressing can be used with compression therapy and is designed to manage varying levels of exudate, from low to high. It is gentle even on fragile skin and does not cause damage to the wound, periwound or surrounding skin. It takes on the toughest challenges impacting adhesion and wear time.

- Advanced adhesive allows for “lift and check” monitoring of wounds without loss of adhesive properties (you can lift the dressing and reposition it multiple times).
- Innovative spoke delivery system enables accurate and easy application in challenging locations like the sacrum.
- Unique multi-layer design absorbs and evaporates moisture while adapting to changing levels of exudate to maintain an optimal wound healing environment.
- Dressing is soft and conformable with a thin, low-profile edge to help minimize the rolling and lifting that can impact adhesion and wear time.
- Adhesive film backing acts as a protective barrier against external bacteria and viruses.
- Waterproof, it can be worn during showering.

Indicated for:

- Wound management
- Pressure ulcer/injury prevention
- Tube/device cushioning or under compression