

TRUGLYDE ABSORB. SUTURE gauge 2/0 circle 3/8 needle 24mm - 70cm - violet



Product Code	22962
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	1781298
CND	H010101020103
EAN/UPC	8023279229622
GMDN	13908

Description

TRUGLYDE - braided & coated polyglycolic acid

- Gauge: 2/0
- Metric: 3.0
- Colour: Violet
- Circle: 3/8
- Needle profile: reverse cutting
- Needle length: 24 mm
- Suture length: 70 cm

Technical details:

- braided and coated polyglycolic acid
- coating of calcium stearate and polycaprolactone
- sterilized by ethylene oxide
- shelf life of 5 years
- wound Support : 28-30 days
- mass Absorption: 70-90 days

Applications:

- for soft tissue approximation during general, gynecological and orthopedic surgeries

Multilanguage manual and box: GB, FR, IT, ES, PT, DE, GR.