



Product Code	29627
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required

EAN/UPC	8023279296273
---------	---------------

Description

SOFTWARE FOR LED VIDEO COLPOSCOPES (29624-5)

HDMI capture card & professional colposcopy software set including: external capture card, software and driver U-Disc, 2 HDMI cables, HDMI connector, dongle for connecting AC-2311 to the desktop computer (or laptop) & monitor at the same time. The software allows to take and collect images for every inspection in order to analyze, diagnose and detect diseases. It can be used in colposcope inspection and dynamic observation of operation treatment.

Features:

- course for displaying, freezing, capturing and deleting image
- taking and collecting 1 to 100 colour images for every inspection
- operating system: Windows 7, Windows 8, Windows 10, Windows 11

175 standard atlas helpful to contrast, analyze and diagnose

Software mark on images

Mark on pathological changes image, add remark, measure pathological change size and magnify partially or totally

Input of clinic symptom

Printing of 19 kinds of reports, RCI evaluation function

Cases information search and editing

Cases be searched by name, case number, age, date and be added, deleted, edited, memorized and reviewed