

TRUGLYDE ABSORB. SUTURE gauge 2/0 circle 1/2
needle 30mm - 90cm - violet



Product Code	22963
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	1781298
CND	H010101020103
EAN/UPC	8023279229639
GMDN	13908

Description

TRUGLYDE - braided & coated polyglycolic acid

- Gauge: 2/0
- Metric: 3.0
- Colour: Violet
- Circle: 1/2
- Needle profile: round body
- Needle length: 30 mm
- Suture length: 90 cm

- Technical details:
- braided and coated polyglycolic acid
 - coating of calcium stearate and polycaprolactone
 - sterilized by ethylene oxide
 - shelf life of 5 years
 - wound Support : 28-30 days
 - mass Absorption: 70-90 days

- Applications:
- for soft tissue approximation during general, gynecological and orthopedic surgeries

Multilanguage manual and box: GB, FR, IT, ES, PT, DE, GR.