



Product Code	59660
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	1984562
CND	M03050202
GTIN	50707387464269

## Description

### 3M™ Soft-Cast™ Casting Tape

- 3M Code: 82102
- Colour: White
- Dimensions: 5 cm x 3.65 m

A fiberglass casting tape that does not set to a completely rigid state, but remains flexible and soft even when it is completely cured.

- the soft, pliable cast edges improve patient comfort
- can be used in combination with rigid fiberglass splints, as soft edges provide additional comfort

#### Suggested Applications:

- sports Medicine: thumb, wrist and ankle sprains
- pediatric Orthopedics: serial casting for club foot treatment
- general Orthopedics: secondary casting, hybrid casting, corsets
- occupational Therapy: removable splints