



Product Code	59701
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	1949415
CND	M03030201
GTIN	50707387788341

## Description

### 3M™ Synthetic Cast Stockinet™

- 3M Code: MS02
- STOCKINET™ 5 cm x 22.8 m

3M™ Stockinet is a synthetic tubular knit made of polyester fibers, used as protection for the patient's skin under synthetic plaster braces.

It is not a sterile product but can be steam sterilized in an autoclave.

Water repellent, making it suitable for use with Soft Cast as a waterproof cast. Suitable for use with all casting materials, including plaster and synthetic.