


PAP TEST KIT - sterile



Product Code	29704
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	I S

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2671276
CND	U0899
EAN/UPC	8023279297041
PARAF	977812940
GMDN	44061

Description

Pap-Test is a painless examination to detect the presence of any anomalous cells on the cervix. It allows an early diagnosis of uterus neck cancer.
The functional configuration of the kit allows a ready to go operation.

The kit includes: one pair latex glove, one side-screw speculum (medium size), one GIMA brush ball, one Ayre spatula A size (wooden), two slides.

Individually packed in sterile pouch.

Multilanguage label and manual: GB, FR, IT, ES, PT, DE, NL, DK, SE, NO, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic.