

INTRAUTERINE DEVICE TCU380 PLUS - sterile



Product Code	29710
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	III

RDM (NSIS)	1509479
CND	U089007
EAN/UPC	3662914002968
PARAF	977812953

Description

INTRAUTERINE DEVICE TCU380 PLUS - sterile

Sterile disposable medical device Class III for intrauterine contraception in women of child bearing age. Can be left inserted for a maximum of 5 years. Packed in a pouch with an insertion tube equipped with a movable flange to aid in gauging the right depth of insertion and a rod. Complies with ISO-7349.

Size: 33x45 mm