

COTTON SWABS - 20 cm



Product Code	29750
Unit of sale	box of 400 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2694991
CND	U0899
EAN/UPC	8023279297508
PARAF	903142026

Description

Cotton Swab

- Shaft material: plastic
- Ø Cotton: 1 cm
- Shaft length: 20 cm

Packed in plastic bags of 25 pcs., in carton of 400 pcs.

Multilingual instructions: GB, FR, IT,ES, PT