

COTTON SWABS - cotton 10 mm (10 plastic bags of 50)



Product Code	29756
Unit of sale	box of 500 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2733131
CND	M9099
EAN/UPC	8023279297560
GMDN	33722

Description

Cotton Swab

- Shaft material: wood
- Ø Cotton: 1 cm
- Shaft length: 15 cm

Packed in plastic bags of 50 pcs., in carton of 500 pcs.