

# MALE URINAL BOTTLE WITH HANDLE 1000 ml



Product Code	25977
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2890244
CND	V0812
GTIN	15060098135420

## Description

### MALE URINAL BOTTLE WITH HANDLE 1000 ml - graduated

Hand-held male urinal bottle made of strong and durable polyethylene.  
Secure fitting attached cap to prevent spills, reduce unpleasant odours and risk of cross infection.  
It can be hand washed or washed in washer disinfectors at 80°C / 90°C; 1000 cycles guaranteed.  
Smooth material surfaces and hollow handle, ensure ease of cleaning.

Size: 180x100x85 mm.  
Easy grip handle

Made in the UK.