

FILE RASP DOUBLE SIDED - 24 cm - metal handle



Product Code	39794
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8023279397949
---------	---------------

PARAF	977816267
-------	-----------

Description

File rasp double sided - 24 cm - metal handle

- STAINLESS STEEL