

RIESTER CUFF 2 TUBES - child



Product Code	49821
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes
EAN/UPC	4045396001077

Description

RIESTER ORIGINAL CUFF - LATEX FREE

- 2 TUBES
- Child velcro cuff
- Riester code: 130
- Colour: light blue
- Size cm: 35.5x10
- Circumference cm: 13 - 17

Compatibility:
All Riester twin tube sphygmomanometers including BigBen