



Product Code	49828
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes
EAN/UPC	4045396168176
PARAF	987475771

## Description

### ONE PIECE RIESTER ORIGINAL CUFF - LATEX FREE

- 2 TUBES - length 50 cm
- Adult large one-piece cuff
- Riester code: 153
- Colour: dark green
- Size cm: 63x16
- Circumference cm: 32 - 41

Compatibility:  
All Riester twin tube sphygmomanometers including BigBen

They prevent cross contaminations. Comfortable soft surface, latex free.  
Washable at 60°C and disinfectable in liquid, reusable up to 10,000 times. ISO 105 E01 tested.