

TRUGLYDE FAST ABSORB. SUTURE gauge 3/0 circle
3/8 needle 20mm - 75cm - undyed



Product Code	22984
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III
UK-REP	No
CH-REP	No
RDM (NSIS)	1920380
CND	H010101020103
EAN/UPC	8023279229844
GMDN	13908

Description

TRUGLYDE FAST - Braided and coated Polyglycolic Acid fast Undyed

- Gauge: 3/0
- Metric: 2.0
- Colour: Undyed
- Circle: 3/8
- Needle profile: reverse cutting
- Needle length: 20 mm
- Suture length: 75 cm

Technical details:

- braided and coated
- coated with calcium stearate and polycaprolactone
- sterilized by ethylene oxide
- shelf life of 5 years
- wound support: 7-10 days
- mass absorption: 40-42 days

Applications:

- for episiotomy
- for circumcision
- closure of skin

Multilanguage manual and box: GB, FR, IT, ES, PT, DE, GR.