



Product Code	49890
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
--------	-----

RDM (NSIS)	2686680
CND	Z1203020501
EAN/UPC	4015672111318

Description

OMRON PROFESSIONAL B.P. MONITOR HBP-1120-E

The HBP-1120 is designed for use in professional settings and is clinically proven to produce fast and reliable results. Irregular heartbeats detection. It is intended for intensive daily use by healthcare professional in hospitals and medical facilities. Zero Indicator function shows that the device is ready for next measurement. It works with AC adaptor or batteries. Class IIa.

Supplied with cuff 22-32 cm, AC adapter, guarantee card and multilingual manual in GB, FR, IT, ES, DE, NL, TR, RU, Arabic.

Technical Specifications

Method: oscillometric, auscultatory
 Measurement range: pressure: 0-300 mmHg; ± 3 mmHg
 pulse rate: 40-200 beats/min.; $\pm 5\%$
 Operating mode: continuous operation
 Display: 7 segment LCD
 Power supply: A C adapter
 Input voltage: A C 100-240 V; 50/60 Hz
 Output voltage: DC 6V ± 0.5 V
 Size - weight: 130x120xh 175 mm - 510 g