

# BLUE SILICONE PESSARY diameter 70 mm - sterile



Product Code	29905
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	No
CH-REP	No

RDM (NSIS)	1295793
CND	U089005
EAN/UPC	8055060761684
PARAF	977813207
GMDN	35237

## Description

Blue silicone pessary, for uterine prolapse.

- Sterile
- Latex free
- Autoclavable at 121°C/15 min or 134°C/10 min for up to 3 times

Individually packed in cardboard box. Multilingual box and instructions: GB, FR, IT, ES, PT, DE, GR.