



| | |
|---------------|----------------|
| Product Code | 49910 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | II A |

| | |
|------------|---------------|
| UK-REP | No |
| CH-REP | No |
| RDM (NSIS) | 2556616 |
| CND | Z1203020501 |
| EAN/UPC | 4015672111745 |
| PARAF | 979256942 |

Description

OMRON M6 COMFORT - HEM-7360-E with Intelli Wrap Cuff 22-42 cm.

- Clinically validated: to verify validations, please see www.dablededucational.org
- Intellisense Technology
- Irregular heartbeat detector: AFib
- Cuff wrap guide
- Intelli Wrap Technology (360° accuracy: accurate measurements from any position around the upper arm)
- Cuff: Intelli wrap cuff 22-42 cm
- Memory measurement: 100 x 2 users + guest mode
- Body movement detection
- Easy high blood pressure colour indicator: signal in case the values exceed the E.H.S. standards for domiciliary auto-measurements (85 systolic pressure and/or 135 diastolic pressure)
- Hypertension
- Average last 3 measurements (Average of 3 last measurements made within 10 min)
- 3 consecutive measurements in AFib mode
- Measurements comparison versus prior readings

Intelli Wrap Cuff fits normal and large arms eliminating inaccurate readings associated with incorrect positioning. AFib indication function.

The Intellisense™ technology ensures that the right inflation pressure is automatically detected.

It detects irregular pulses and automatically decides if the measurement result is reliable or needs to be repeated.

Only accurate results are displayed.

Multilingual box and manual: GB, FR, IT, ES, DE, NL, RU, TR, Arabic.

Technical Specifications

Measuring method: oscillometric
 Measuring range:
 - Pressure: 0 mmHg to 299 mmHg \pm 3 mmHg
 - Pulse: 40 to 180 min \pm 5% of display reading
 Memory (with date and time): 100x2 users
 Inflation: fuzzy logic controlled by electric pump
 Deflation: automatic pressure release valve
 Rating: 6VCC 4W
 Power source: 4 "AA" batteries 1.5 V
 Battery life (N° of measurements): 1000

Protection against electric shock: internally powered ME equipment (when using batteries)

Weight (console+cuff) without batteries: 460+163 g

Size: 191x120xh 85 mm

Cuff materials: nylon and polyester

Measurable circumference (cm): 22-42 (arm)