



Product Code	49911
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2576435
CND	Z1203020501
EAN/UPC	4015672113978
PARAF	988050086

## Description

### OMRON M6 COMFORT - AFIB - HEM-7380-E

With Intelli Wrap Cuff 22-42 cm.

It detects the possibility of Atrial fibrillation, enabling home monitoring. It also takes 3 consecutive readings at 30 second intervals, then displays the average to give more accurate indication of blood pressure.

The Intellisense™ technology ensures that the right inflation pressure is automatically detected. It detects irregular pulses and automatically decides if the measurement result is reliable or needs to be repeated. Only accurate results are displayed.

Clinically validated (To verify validations, please see [www.dablededucational.org](http://www.dablededucational.org))

Intellisense technology

Irregular heartbeats detector (AFib)

Cuff wrap guide

Intelli Wrap Technology (360° accuracy: accurate measurements from any position around the upper arm)

Pre-formed cuff: model Intelli wrap / size (cm) 22-42 (Pre-formed cuff, excluding M2+, M2 Intelli IT, M3)

Memory measur. x user 100x2

Guest mode

Body movement detection

Easy high blood pressure colour indicator (Signal in case the values exceed the E.H.S. standards for domiciliary auto-measurements (85 systolic pressure and/or 135 diastolic pressure))

Hypertension

Average last 3 measurements (Average of 3 last measurements made within 10 min)

AFib detection in each measurement

Measurements comparison versus prior readings

## Technical Specifications

Measuring method: oscillometric

Measuring range:

- Pressure: 0 mmHg to 299 mmHg  $\pm$  3 mmHg

- Pulse: 40 to 180 min  $\pm$  5% of display reading

Memory (with date and time): 100x2 users

Inflation: fuzzy logic controlled by electric pump

Deflation: automatic pressure release valve

Rating: 6VCC 4W

Power source: 4 "AA " batteries 1.5 V

Battery life (N° of measurements): 1000

Protection against electric shock: internally powered ME equipment (when using batteries)

Weight (console+cuff) without batteries: 460+163 g

Size: 191x120xh 85 mm

Cuff materials: nylon and polyester

Measurable circumference (cm): 22-42 (arm)