

EVOLUTION DOUBLE HYSTEROMETER



Product Code	29929
Unit of sale	box of 25 pcs
Minimum order	1
Type	Medical device
Class	I S

RDM (NSIS)	2468031
CND	U089003
EAN/UPC	8023279299298

Description

EVOLUTION DOUBLE HYSTEROMETER 30 cm

The most versatile, economical and eco-friendly single-use hysterometer.
Estimating the uterine depth, checking cervical canal permeability and womb orientation.
Sterile, disposable medical device, class I.
Latex free, PVC free.
Shelf life: 5 years.
Multilingual manual: GB, FR, IT, ES, PT, DE.

Innovative double track system.
2 flexible and shape memory tips with gradual dilation.
For all types of IUD insertion and womb shapes.

- One side:
Diameter: from CH 10 to CH 14
(from 3.3 mm to 4.7 mm)

- One side:
Diameter: from CH 8 to CH 12
(from 2.66 mm to 4 mm)

Pale pink color allowing better observation of cervical mucus.