

FARMAGEL C LUBRICATING GEL WITH
CHLOREXEDINA 12.5 g - sterile



Product Code	29961
Unit of sale	box of 25 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	1420502
CND	V9007
EAN/UPC	8004007001673
GMDN	37717

Description

FARMAGEL® C LUBRICATING GEL WITH CHLORHEXIDINE - single use squeeze bulb 12.5 g - sterile - box of 25

Farmagel® C is a sterile water-soluble, grease-free, transparent lubricating gel.
Suitable for endoscopic procedures which require lubricating action.
It reduces friction between probe and mucosa facilitating the entry of diagnostic or therapeutic devices into the body.
It does not irritate and is easily removable from skin and clothing.

Ingredients: chlorhexidine (0.05%), hydroxyethyl cellulose, glycerin, purified water, hydrochloric acid. Class IIa medical device.

Made in Italy.