



GIMA

1 pz.
1 pc.



GIMA 23456

LOT 000005

RILEVATORE DI VENE PROFESSIONALE QV-500
QV-500 PROFESSIONAL VEIN FINDER
DÉTECTEUR DE VEINES PROFESSIONNEL QV-500
DETECTOR DE VENAS PROFESIONAL QV-500



SN

UDI



Qualmedi Technology Co., Ltd.

A302 Room, No.23, HangBu Road, Feixi County Economic
Development Zone, Hefei City, Anhui Province, China.

sales@qualmedi.com

Made in China

REF

QV-500

EC

REP

Kingsmead Service B.V.
Zonnehof 36, 2632 BE, Nootdorp,
Netherlands



Gima S.p.A.

Via Marconi, 1 – 20060 Gessate (MI) – Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com



MD



Rev.1.10.24

