



GIMA

1 pz. Rev.2.03.25
1 pc.



GIMA 28580

LOT 000003

KIT 8010 MATERASSO A PRESSIONE ALTERNATA NYLON PVC+POMPA
KIT 8010 NYLON PVC ALTERNATING PRESSURE MATTRESS + PUMP
KIT 8010 MATELAS À PRESSION ALTERNÉE NYLON PVC + POMPE
KIT 8010 COLCHÓN PRESIÓN ALTERNADA NYLON PVC + BOMBA



Guangdong Yuehua Medical Instrument Factory Co., Ltd.
Rongsheng Science and Technology Zone, Daxue Road, Shantou,
Guangdong, China
Made in China

REF QDC-8010 + P3000A2QB



UDI

EC REP

Eunitor GmbH

Kennedydamm, 5, Düsseldorf, 40476, Germany

Importato da/Imported by/Importé par/Importado por
Gima S.p.A
Via Marconi, 1 – 20060 Gessate (Mi) – Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

