



GIMA

1 pz. Rev.2.03.25
1 pc.

GIMA 28582 **LOT** 000003

KIT 8080 MATERASSO A PRESSIONE ALTERNATA NYLON TPU+POMPA
KIT 8080 NYLON TPU ALTERNATING PRESSURE MATTRESS + PUMP
KIT 8080 MATELAS À PRESSION ALTERNÉE NYLON TPU + POMPE
KIT 8080 COLCHÓN PRESIÓN ALTERNADA NYLON TPU + BOMBA



Guangdong Yuehua Medical Instrument Factory Co., Ltd.
Rongsheng Science and Technology Zone, Daxue Road, Shantou,
Guangdong, China
Made in China

REF QDC-8080+P3000A2QB3




Eunitor GmbH

Kennedydamm, 5, Düsseldorf, 40476, Germany

Importato da/Imported by/Importé par/Importado por
Gima S.p.A

Via Marconi, 1 – 20060 Gessate (Mi) – Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com



8

023279

285826



MD



UDI