

REF 830-340-30 (GIMA 30826)

LOT 000019



1 pz.



CAVO A FIBRE OTT. 4,8x3.000 mm - senza adattatori
F.O.CABLE 4.8x3000 mm - without adapters
CÂBLE EN F.O. 4,8x3000 mm - sans adaptateur
CABLE DE FIBRA OPT. 4,8x3.000 mm - sin adaptado



Medical Devices (Pvt) Ltd
Wazirabad Road, Ugoki Sialkot. Pakistan
Made in Pakistan



Obelis s.a. Bd Général Wahis 53
B-1030 Brussels, Belgium



GIMA S.p.A. - Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

