



1 pz.



REF 31427

LOT 000016

LAMPADINA OTOSCOPIO XENON-ALOGENA - ricambio
XENON-HALOGEN OTOSCOPE BULB - spare
AMPOULE XÉNON HALOGÈNE POUR OTOSCOPE - rechange
BOMBILLA otoscopio XENON-HALÓGENA - de recambio



Importato da/Imported by/Importé par/Importado por:
GIMA S.p.A. - Via Marconi, 1 - 20060 Gessate (MI) - Italy
Made in Pakistan - www.gimaitaly.com