

GIMA

Conf 250 pz.
Box of 250 pcs.

REF 31488

LOT 000040

MINI SPECULUM MONOUSO diam. 4 mm - neri
MINI EAR SPECULUM diam. 4 mm - black
MINI-SPÉCULUM AURICULAIRE diam. 4 mm - noir
MINI ESPÉCULO diám. 4 mm - negro - en dispensador



8 023279 314885



PAP



PAP



Rev. 3-11/2024



Gima S.p.A. - Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com - Made in China

